

| TITMUS VSERIES MODEL COMPARISON   |               |             |               |
|-----------------------------------|---------------|-------------|---------------|
| Model                             | Titmus V1     | Titmus V2   | Titmus V4     |
| <b>Control Panel</b>              | ✓             |             | ✓             |
| <b>Number of Slides</b>           | 4             | 8           | 8             |
| <b>No Moving Parts</b>            | ✓             |             |               |
| <b>Ergonomics</b>                 | ✓             | ✓           | ✓             |
| <b>Light Source</b>               | LED           | Fluorescent | Fluorescent   |
| <b>Portability</b>                | ✓             | ✓           | ✓             |
| <b>Soft Foam Headrest</b>         | ✓             | ✓           | ✓             |
| <b>Head Sensor</b>                | ✓             |             | ✓             |
| <b>Peripheral Vision Testing</b>  | ✓             | ✓           | ✓             |
| <b>Test Slide Advance/Reverse</b> | Control Panel | Manual      | Control Panel |
| <b>Night Vision Testing</b>       | ✓             |             | ✓             |
| <b>Glare Recovery Testing</b>     | ✓             |             |               |
| <b>Carrying Case</b>              | Optional      | Optional    | Optional      |
| <b>Standard Accessories</b>       | ✓             | ✓           | ✓             |

#### HOW TO SELECT A VISION SCREENER

**STEP 1:** Select Model Titmus V1 or Titmus V4 based on features & benefits, ex. Titmus V1

**STEP 2:** Select a slide set based on the type of testing, ex. DMV, Professional, etc.

**STEP 3:** Provide the Vision Screener Model AND Slide set information to your distributor, ex. Titmus V1 DMV